

Product Datasheet

Human SNRPC protein orb428518

Description	Recombinant of human SNRPC protein
Reactivity	Human
Conjugation	Unconjugated
Target	SNRPC
Form/Appearance	Sterile Filtered clear solution. SNRPC protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 50% glycerol, 0.3M NaCl, 5mM DTT and 2mM EDTA.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMPKFYCD YCDTYLTHDS PSVRKTHCSG RKHKENVKDY YQKWMEEQAQ SLIDKTTAAF QQGKIPPTPF SAPPAGAMI PPPPSLPGPP RPGMMPAPHM GGPPMPMMG PPPPGMPVG PAPGMRPPMG GHMPMPGPP MMRPPARPM VPTRPGMTRP DR.
Purity	>85% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.