


Product Datasheet

Human Pro-MMP13 ELISA Kit orb397082

Description	Sandwich ELISA for quantitative detection of Human Pro-MMP13 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	30 pg/ml
Tested Applications	ELISA
Target	MMP13
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma



Representative Standard Curve