

## Product Datasheet

### Human SPECC1 protein orb428566

<b>Description</b>	Recombinant of human SPECC1 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	SPECC1
<b>Form/Appearance</b>	Sterile Filtered colorless solution. SPECC1 protein solution (0.25mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 10% glycerol 1mM DTT and 0.1M NaCl.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSTQEGKII ELEQKCTGIL EQGRFEREKL LNIQQQLTCS LRKVEEENQG ALEMIKRLKE ENEKLNFEFLERHNNNMMA KTLIEECRVTL EGLKMENGSL KSHLQGEKQK ATEASAVEQT AESCEVQEML KVARAEKDLL ELSCNELRQE LLKANGEIKH VSSLLAKVEK DYSYLKEICD HQAEQLSRTS LKLQEKASES DAEIKDMKET IFELEDQVEQ HRAVKLHNNQ LISELESSVI KLEEQKSDLE RQLKTLTKQM KEETEEWRRF QADLQTAVVV ANDIKCEAQQ ELRTVTKRLL EEEEKNA.
<b>Purity</b>	>80.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.