

Product Datasheet

Human SUB1 protein orb428665

Description	Recombinant of human SUB1 protein
Reactivity	Human
Conjugation	Unconjugated
Target	SUB1
Form/Appearance	Sterile filtered colorless solution. SUB1 Human solution containing 20mM Tris-HCl pH-8, 0.2M NaCl, 5mM DTT & 20% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MPKSKELVSS SSSGSDSDSE VDKLLKRRKKQ VAPEKPVKKQ KTGETSRALS SSKQSSSSRD DNMFQIGKMR YVSVRDFKGG VLIDIREYWM DPEGEMKPGR KGISLNPEQW SQLKEQISDI DDAVRKL.
Purity	>90% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.