

Product Datasheet

T.pallidum p47 (Partial) protein orb428710

Description	Recombinant of T.pallidum p47 (Partial) protein
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	T.pallidum p47 (Partial)
Form/Appearance	70mM Tris-HCl pH8.0, 50mM NaCl, 50% Glycerol and 1.5M urea.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Purity	>95% pure as determined by 10% PAGE (coomassie staining).
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.