

Product Datasheet

Human TARBP2 protein orb428734

Description	Recombinant of human TARBP2 protein
Reactivity	Human
Conjugation	Unconjugated
Target	TARBP2
Form/Appearance	Sterile Filtered colorless solution. TARBP2 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 10% glycerol and 0.4M Urea.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MLAANPGKTP ISLLQEYGTR IGKTPVYDLL KAEGQAHQPN FTFRVTVGDT SCTGQGPSKK AAKHKAAEVA LKHLKGG SML EPALEDSSSF SPLDSSLPED IPVFTA AAAA TPVPSVVLTR SPPMELQPPV SPQQSECNPV GALQELVVQK GWRLPEYTVT QESGPAHRKE FTMTCRVERF IEIGSGTSKK LAKRNAAAKM LLRVHTVPLD ARDGNEVEPD DDHFSIGVGS RLDGLRNRGP GCTWDSL RNS VGEKILSLRS CSLGSLGALG PACCRVLSSEL SEEQAFHVS Y LDIEELSLSG LCQCLVELST QPATVCHGSA TTREAARGEA ARRALQYLKI MAGSK.
Purity	>85.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.