

Product Datasheet

Human TIE1 Fc protein orb428846

Description	Recombinant of human TIE1 Fc protein
Reactivity	Human
Conjugation	Unconjugated
Target	TIE1 Fc
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder. TIE-1 Fc Chimera was lyophilized from a concentrated (1mg/ml) sterile solution containing 20mM Tris, 0.5M NaCl, 10% Sucrose.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Purity	>90.0% as determined by:(A) Analysis by RP-HPLC.(B) Analysis by SDS-PAGE.
Source	Insect Cells.
Solubility (25°C)	It is recommended to reconstitute the lyophilized TIE-1 Fc Chimera in sterile water not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.