


Product Datasheet

Mouse HIF1 alpha ELISA Kit orb397111

Description	Sandwich ELISA for quantitative detection of Mouse HIF1 alpha concentrations.
Reactivity	Mouse
Range	125 pg/ml - 8000 pg/ml
Sensitivity	60 pg/ml
Tested Applications	ELISA
Target	HIF1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma, Tissue Homogenates



Representative Standard Curve