

## Product Datasheet

### Human TNNI3 (His) protein orb428931

<b>Description</b>	Recombinant of human TNNI3 (His) protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	TNNI3 (His)
<b>Form/Appearance</b>	Sterile Filtered clear solution. TNNI3 protein solution in PBS and 25mM arginine.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	>95% pure as determined by 12% PAGE (coomassie staining).
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.