

Product Datasheet

Human TRPT1 protein orb429004

Description	Recombinant of human TRPT1 protein
Reactivity	Human
Conjugation	Unconjugated
Target	TRPT1
Form/Appearance	Sterile Filtered colorless solution. The TRPT1 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SGLVPRGSH MGSHMNFSGG GRQEAAGSRG RRAPRPREQD RDVQLSKALS YALRHGALKL GLPMGADGFV PLGTLQLLPQ FRGFAEDVQ RVVDTNRKQR FALQLGDPST GLLIRANQGH SLQVPKLELM PLETPQALPP MLVHGTFWKH WPSILLKGLS CQGRTHIHLA PGLPGDPGII SGMRSHEIA VFIDGPLALA DGIPFFRSAN GVILTPGNTD GFLLPKYFKE ALQLRPRTRKP LSLAGDEETE CQSSPKHSSR ERRRIQQ.
Purity	>95.0% as determined by SDS-PAGE.
Source	Escherichia Coli.
Expiration Date	6 months from date of receipt.