

## Product Datasheet

### Human UBA5 protein orb429061

<b>Description</b>	Recombinant of human UBA5 protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	UBA5
<b>Form/Appearance</b>	Sterile filtered colorless solution. The UBA5 solution (0.5mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 50mM NaCl.
<b>Storage</b>	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	MGSSHHHHHH SSGLVPRGSH MGSHEMAESVE RLQQRVQELE RELAQERSLQ VPRSGDGGGG RVRIEKMSSE VVDSNPYSRL MALKRMGIVS DYKIRTFV AIVGVGGVGS VTAEMLTRCG IGKLLLFDYD KVELANMNRL FFQPHQAGLS KVQAAEHLR NINPDVLFV HNYNITTVENFQHFMDRISN GGLEEGKPV LVLSCVDNFE ARMTINTACN ELGQTWMESE VSENAVSGHI QLIIPGESAC FACAPPLVVA ANIDEKTLKR EGVCAASLPT TMGVVAGILV QNVLKFLNFTVSVFYLGN AMQDFFPTMS MKPNPQCDDR NCRKQEEYK KKVAALPKQE VIQEEEEIIH EDNEWGIELV SEVSEELKN FSGPVPDLPE GITVAYTIPK KQEDSVTELT VEDSGESLED LMAKMKNM.
<b>Purity</b>	>90.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli.
<b>Expiration Date</b>	6 months from date of receipt.