


Product Datasheet

Mouse MMP8 ELISA Kit orb397127

Description	Sandwich ELISA for quantitative detection of Mouse MMP8 concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Sensitivity	80 pg/ml
Tested Applications	ELISA
Target	MMP8
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates



Representative Standard Curve