

Product Datasheet

Human VAMP7 protein orb429164

Description	Recombinant of human VAMP7 protein
Reactivity	Human
Conjugation	Unconjugated
Target	VAMP7
Form/Appearance	Sterile Filtered colorless solution. The VAMP7 solution (0.5mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl and 20% glycerol.
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only.
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MGSMALFAV VARGTTILAK HAWCGGNFLE VTEQILAKIP SENNKLTYSYH GNYLFHYICQ DRIVYLCITD DDFERSRAFN FLNEIKKRFQ TTYGSRAQTA LPYAMNSEFS SVLAAQLKHH SENKGLDKVM ETQAQVDELK GIMVRNIDL AQRGERLELL IDKTENLVDS SVTFKTTSRN LARAMCMKNL K
Purity	>90% as determined by SDS-PAGE.
Source	E.coli.
Expiration Date	6 months from date of receipt.