


## Product Datasheet

### Mouse TROY ELISA Kit orb397133

<b>Description</b>	Sandwich ELISA for quantitative detection of Mouse TROY concentrations.
<b>Reactivity</b>	Mouse
<b>Range</b>	62.5 pg/ml - 4000 pg/ml
<b>Sensitivity</b>	30 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	TNFRSF19
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<b>Detection Method:</b> Colorimetric <b>Detection Range:</b> Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates



Representative Standard Curve