


Product Datasheet

Rat CXCL3 ELISA Kit orb397137

Description	Sandwich ELISA for quantitative detection of Rat CXCL3 concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Sensitivity	15 pg/ml
Tested Applications	ELISA
Target	CXCL3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma



Representative Standard Curve