


Product Datasheet

Human MMP10 ELISA Kit orb402645

Description	Sandwich High Sensitivity ELISA kit for Quantitative Detection of Human MMP10 ELISA Kit
Reactivity	Human
Range	78pg/ml-5000pg/ml
Sensitivity	
Tested Applications	ELISA
Target	MMP10
Storage	Store at 4°C for 6 months and at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice).
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<p>Typical Data Sample Data: Concentration(pg/ml) 0 78 156 313 625 1250 2500 5000 O.D. 0.038 0.147 0.215 0.349 0.564 0.961 1.66 2.124 TMB Incubation Time: 15-20min</p>



Human MMP-10 PicoKine
ELISA
Kit standard curve