


Product Datasheet

Mouse Ctf1 ELISA Kit (orb402667)

Description	Sandwich High Sensitivity ELISA kit for Quantitative Detection of Mous Ctf1 ELISA Kit...
Reactivity	Mouse
Range	156pg/ml-10000pg/ml
Sensitivity	< 10pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S2-A203
Target	Ctf1
Storage	Store at 4°C for 6 months and at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice).
Note	For research use only.
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q60753
Expiration Date	Please enquire.
Application Notes	<p>Typical Data</p> <p>Sample Data: Concentration(pg/ml) 0 156 312 625 1250 2500 5000 10000 O.D. 0.006 0.117 0.212 0.381 0.637 1.064 1.476 1.917</p> <p>TMB Incubation Time: 20-25min</p>



Mouse Cardiotrophin-1 PicoKine ELISA Kit standard curve