


Product Datasheet

Human DAG1 ELISA Kit orb402728

Description	Sandwich High Sensitivity ELISA kit for Quantitative Detection of Human DAG1 ELISA Kit
Reactivity	Human
Range	156pg/ml-10000pg/ml
Sensitivity	
Tested Applications	ELISA
Target	DAG1
Storage	Store at 4°C for 6 months and at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice).
Note	For research use only.
Uniprot ID	Q14118
Entrez	1605
Expiration Date	Please enquire.
Application Notes	Typical Data Sample Data: Concentration(pg/ml) 0 156 312 625 1250 2500 5000 10000 O.D. 0.008 0.069 0.122 0.229 0.452 0.812 1.498 2.120 TMB Incubation Time: 20-25min



Human Dystroglycan
PicoKine
ELISA Kit standard curve