


Product Datasheet

SARS-CoV-2 trimeric soluble full-length Spike protein, B1.617 cluster, S protein orb864239

Description	Recombinant human SARS-CoV-2 trimeric soluble full-length Spike protein, B1.617 cluster, S protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Immunogen	Unconjugated
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	8xHis at the C-terminus
Note	For research use only.
Protein Sequence	16-1211
Purity	> 95%
MW	137
Source	HEK293 cells
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.



SARS-CoV-2 trimeric soluble full-length Spike protein, B1.617 cluster, S protein