

## Product Datasheet

# Viral HEPATITIS B CORE ANTIGEN protein (orb435682)

<b>Description</b>	Recombinant of viral HEPATITIS B CORE ANTIGEN protein...
<b>Reactivity</b>	Virus
<b>Tested Applications</b>	ELISA
<b>Target</b>	HEPATITIS B CORE ANTIGEN
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only.
<b>Clonality</b>	Recombinant
<b>Purity</b>	Rec. Protein
<b>Expiration Date</b>	6 months from date of receipt.