

Product Datasheet

Bovine AHSG ELISA Kit orb403036

Description	Bovine Alpha-2-HS-glycoprotein(AHSG) ELISA kit
Reactivity	Bovine
Range	125 ng/ml-8000 ng/ml
Sensitivity	31.25 ng/ml
Tested Applications	ELISA
Target	AHSG
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma , tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm