

Product Datasheet

Bovine ANGPTL4 ELISA Kit [Out of Stock] orb403038

Description	Bovine Angiopoietin-related protein 4(ANGPTL4) ELISA kit
Reactivity	Bovine
Sensitivity	0.027 ng/ml
Tested Applications	ELISA
Target	ANGPTL4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm