

Product Datasheet

Bovine MSTN ELISA Kit orb403086

Description	Bovine Growth/differentiation factor 8(MSTN) ELISA kit
Reactivity	Bovine
Range	6.25 ng/mL-400 ng/mL
Sensitivity	1.56 ng/mL
Tested Applications	ELISA
Target	MSTN
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O18836
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm