

Product Datasheet

Bovine IFNW1 ELISA Kit orb403098

| | |
|----------------------------|---|
| Description | Bovine Interferon omega-1(IFNW1) ELISA kit |
| Reactivity | Bovine |
| Sensitivity | 15.6 pg/ml. |
| Tested Applications | ELISA |
| Target | IFNW1 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |