

Product Datasheet

Bovine IL2RA ELISA Kit [Out of Stock] orb403107

Description	Bovine Interleukin-2 receptor subunit alpha(IL2RA) ELISA kit
Reactivity	Bovine
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	IL2RA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm