

Product Datasheet

Bovine MB ELISA Kit orb403126

Description	Bovine Myoglobin(MB) ELISA kit
Reactivity	Bovine
Range	15.6 ng/ml-1000 ng/ml
Sensitivity	3.9 ng/ml
Tested Applications	ELISA
Target	MB
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue Homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm