

Product Datasheet

Bovine TYR ELISA Kit [Out of Stock] orb403172

Description	Bovine Tyrosinase(TYR) ELISA kit
Reactivity	Bovine
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.8 pg/ml
Tested Applications	ELISA
Target	TYR
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm