

Product Datasheet

Bovine VEGFA ELISA Kit (orb403174)

Description	Bovine Vascular endothelial growth factor A(VEGFA) ELISA kit
Reactivity	Bovine
Range	47 pg/mL-3000 pg/mL
Sensitivity	11.75 pg/mL
Tested Applications	ELISA
Target	VEGFA
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P15691
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm