

## Product Datasheet

### 10x Phosphate Buffered Saline orb436759

<b>Description</b>	10x Phosphate Buffered Saline
<b>Hazard Information</b>	Non-Toxic
<b>Target</b>	10x Phosphate Buffered Saline
<b>Form/Appearance</b>	Concentrate
<b>Storage</b>	Store at 4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Note</b>	For research use only.
<b>Purity</b>	10 x concentrate pH7.3 +/-0.1
<b>Expiration Date</b>	12 months from date of receipt.