

## Product Datasheet

### Bovine DEFB1 ELISA Kit orb403194

|                            |                                                                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Bovine Beta-defensin 1(DEFB1) ELISA kit                                                                                                                                                                                           |
| <b>Reactivity</b>          | Bovine                                                                                                                                                                                                                            |
| <b>Range</b>               | 3.12 ng/ml-200 ng/ml                                                                                                                                                                                                              |
| <b>Sensitivity</b>         | 0.78ng/ml                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                             |
| <b>Target</b>              | DEFB1                                                                                                                                                                                                                             |
| <b>Note</b>                | For research use only.                                                                                                                                                                                                            |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                                   |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |