

Product Datasheet

Bovine CSN2 ELISA Kit orb403196

Description	Bovine Beta-casein(CSN2) ELISA kit
Reactivity	Bovine
Range	14 µg/ml-900 µg/ml
Sensitivity	7 µg/ml
Tested Applications	ELISA
Target	CSN2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: breast milk• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm