

Product Datasheet

Bovine β -OHB ELISA Kit orb403198

Description	Bovine Beta-Hydroxybutyric Acid, β -OHB ELISA Kit
Reactivity	Bovine
Range	50 nmol/ml-800 nmol/ml
Sensitivity	9.5 nmol/ml
Tested Applications	ELISA
Target	β -OHB
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm