

Product Datasheet

Bovine IL6 ELISA Kit orb403210

Description	Bovine Interleukin 6,IL-6 ELISA Kit
Reactivity	Bovine
Range	5 pg/ml-1000 pg/ml
Sensitivity	2.5 pg/ml
Tested Applications	ELISA
Target	IL6
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm