

Product Datasheet

Bovine FGF21 ELISA Kit orb403219

Description	Bovine fibroblast growth factor 21 (FGF21) ELISA kit
Reactivity	Bovine
Range	62.5 pg/ml-4000 pg/ml
Sensitivity	15.6 pg/ml
Tested Applications	ELISA
Target	FGF21
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm