
Product Datasheet

Bovine FGF21 ELISA Kit (orb403219)

Description	Bovine fibroblast growth factor 21 (FGF21) ELISA kit
Reactivity	Bovine
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Tested Applications	ELISA
Target	FGF21
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm