

Product Datasheet

Bovine ACTH ELISA Kit orb403223

Description	Bovine adrenocorticotrop hormone (ACTH) ELISA kit
Reactivity	Bovine
Range	0.12 pg/ml -500 pg/ml
Sensitivity	0.12 pg/ml
Tested Applications	ELISA
Target	ACTH
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm