

Product Datasheet

Bovine Casp7 ELISA Kit (orb403243)

Description	Bovine Caspase 7(Casp-7) ELISA Kit
Reactivity	Bovine
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Tested Applications	ELISA
Target	Casp7
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	F1MD58
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm