

Product Datasheet

Bovine Caspase 8 ELISA Kit orb403244

Description	Bovine Caspase 8(Casp-8) ELISA Kit
Reactivity	Bovine
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156ng/ml
Tested Applications	ELISA
Target	Caspase 8
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm