

Product Datasheet

Bovine LH ELISA Kit orb403251

Description	Bovine luteinizing hormone(LH)ELISA kit
Reactivity	Bovine
Range	1.25 mIU/ml-100 mIU/ml
Sensitivity	1.56 mIU/ml
Tested Applications	ELISA
Target	LH
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm