

Product Datasheet

Bovine FGF2 ELISA Kit orb403255

Description	Bovine Heparin-binding growth factor 2(FGF2) ELISA kit
Reactivity	Bovine
Range	6.25 pg/ml-400 pg/ml
Sensitivity	1.56 pg/ml
Tested Applications	ELISA
Target	FGF2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm