

Product Datasheet

Bovine Calcitonin ELISA Kit orb403256

Description	Bovine Calcitonin(CALCA/CALC) ELISA kit
Reactivity	Bovine
Range	6.25 pg/ml-400 pg/ml
Sensitivity	25 pg/ml
Tested Applications	ELISA
Target	Calcitonin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm