

Product Datasheet

Bovine IgE ELISA Kit orb403267

Description	Bovine IgE ELISA Kit
Reactivity	Bovine
Range	62.5 ng/mL-4000 ng/mL
Sensitivity	15.6 ng/mL
Tested Applications	ELISA
Target	IgE
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm