

Product Datasheet

Bovine HSP-70 ELISA Kit orb403272

Description	Bovine Heat Shock Protein 70,HSP-70 ELISA Kit
Reactivity	Bovine
Range	1.25 ng/ml-20 ng/ml
Sensitivity	1.25 ng/ml
Tested Applications	ELISA
Target	HSP-70
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm