

Product Datasheet

Bovine HSPA2 ELISA Kit orb403273

Description	Bovine Heat Shock Protein 72 (HSP-72) ELISA Kit
Reactivity	Bovine
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	HSPA2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P34933
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm