

Product Datasheet

Bovine PAG1 ELISA Kit orb403274

Description	Bovine Pregnancy Specific Protein B(PSPB)ELISA Kit
Reactivity	Bovine
Range	0.06 ng/ml-24 ng/ml
Sensitivity	0.06 ng/ml
Tested Applications	ELISA
Target	PAG1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm