

Product Datasheet

Bovine PAG1 ELISA Kit orb403274

Description	Bovine Pregnancy Specific Protein B(PSPB)ELISA Kit
Reactivity	Bovine
Range	0.06 ng/mL-24 ng/mL
Sensitivity	0.06 ng/mL
Tested Applications	ELISA
Target	PAG1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q29432
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm