

Product Datasheet

Bovine LPO ELISA Kit orb403277

Description	Bovine lactoperoxidase (LPO) ELISA kit
Reactivity	Bovine
Range	1.56 mU/ml-100 mU/ml
Sensitivity	1.56 mU/ml
Tested Applications	ELISA
Target	LPO
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: breast milk• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm