
Product Datasheet

Bovine LEP ELISA Kit (orb403284)

Description	Bovine Leptin,LEP ELISA Kit
Reactivity	Bovine
Range	3.12 ng/mL-50 ng/mL
Sensitivity	3.12 ng/mL
Tested Applications	ELISA
Target	LEP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P50595
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm