

Product Datasheet

Bovine FN1 ELISA Kit orb403295

Description	Bovine Fibronectin, FN ELISA Kit
Reactivity	Bovine
Range	0.137 µg/ml -100 µg/ml
Sensitivity	0.137 µg/ml
Tested Applications	ELISA
Target	FN1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm