

Product Datasheet

Bovine IGF1 ELISA Kit orb403306

Description	Bovine Insulin-like growth factor 1,IGF-1 ELISA Kit
Reactivity	Bovine
Range	0.75 ng/mL-60 ng/mL
Sensitivity	0.94 ng/mL
Tested Applications	ELISA
Target	IGF1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P07455
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm